Vlectro Alcohol tester VL-AT

# User Manual

Vlectro VL-AT

## ■ Introduction

Thanks for purchasing our Vlectro VL-AT professional alcohol tester, please study this manual carefully before using.

The VL-AT personal breathalyzer, is a hand held professional personal use semiconductor type alcohol analyzer with LCD digital display.


## ■ Humanized safe high-tech design

◆Key parameters such as sampling analysis, test accuracy and display range designed according to the EN 16280 standard, which drafted and officially adopted by the European Union in the field of breath alcohol testing devices used by the general public in the year of 2012.

◆Single button operation design, easy use.

**Note: Manufacturer reserves the right to modify the contents of this manual without prior notice.**

## ■ Parts Instruction

No.	Item	Function
1	Mouthpiece	Breath test inlet
2	LCD	Display
3	ON/OFF KEY	⏻ Key: ON/OFF and Test
4	Vent outlet	Vent outlet
5	Buzzer	Sound prompt
6	Battery compartment	To install 2 * AAA alkaline batteries

## ■ Operation Instruction

Read the instructions carefully before use.

### 1. How to use

Rinse your mouth with water before using.

Avoid smoking or drinking alcohol for at least 20 minutes before using.

### 2. Insert the batteries

- Open the battery cover on the back of the device.


Vlectro VL-AT

- Insert 2 x 1.5V alkaline batteries (AAA) into the battery compartment.

### 3. Power ON/OFF

#### 3.1 Power ON

After insert batteries, keep press  key for 2,5s until device with “di-di” voice warning, then device enter into the warm-up mode, interface show as following picture

Power on interface

#### 3.2 Power OFF

Under warm-up and Heating state, keep press  for 2,5s, the device will power off;

When the result is displayed, without any key operation within 60s, device will power off automatically.


Power Off interface

### 4. Formal Test

#### 4.1 Intelligent self-checking


Turn on the device, will enter into the warm-up state automatically, to reset the sensor into a better working state, interface as following picture

Warm-up mode

## 4.2 Wait for preheating

After warm-up, will enter into heating mode to active the warming of sensor, countdown 19s (19-00), to make sensor work better and more sensitive.

Countdown interface

### ➤ Tips: Alarm level setting

Under self-checking and preheating state, double press  button to enter alarm level setting mode. Short press the  button again to select alarm level as b2,b3,b5,b8, after setted, waiting for 2 seconds, device will enter into test mode automatically.

b2: Means alarm level as 0.20%BAC(=0.020%BAC or 0.100mg/L)

b3: Means alarm level as 0.30%BAC(=0.030%BAC or 0.150mg/L)


b5: Means alarm level as 0.50%BAC(=0.050%BAC or 0.250mg/L)

b8: Means alarm level as 0.80%BAC(=0.080%BAC or 0.400mg/L)

## 4.3 Blow test


After 19s of heating, when the “blow” flashes, the instrument can be breath tested by the mouthpiece. Press the ON/OFF button one time when starting the breath test, the instrument is officially sampled. “bLo” flashes for about 25 seconds ,Please carry out the measurement within this period.

Blow continuously through the mouthpiece for about 5 seconds, a long beep tone sounds continuously while blowing. Only stop blowing until the beep tone stops together with dynamic progress bar full displayed. Two short “di-di” sounds will be heard when the sample is being taken successfully, test result will be displayed accordingly.

Ready-blow interface

Dynamic progress bar display on LCD while blowing, as following pictures

Dynamic progress bar

Device will power off automatically, when user fails to start the blowing test within 25s.

Test reading interface

◆ Under test reading interface, can press  key to re-test, or device will power off automatically if there's no any operation within 60 seconds.

◆ Under test reading interface, keep press  key to switch units between mg/L,%,%, next test reading will show with this re-setted unit

#### 4.4 Test result show

If the test result exceeds the alert level, the alarm indicator light turns on red and flashes 10 times together with “di-di” sound prompt.

While test result lower than alarm level, the alert indicator light turns on green.

The reading range is 0.00%BAC -2.00%BAC.

“2.00%BAC” will be shown if the reading is higher than 2.00%BAC.

“0.00‰BAC” will be shown if the reading is lower than 0.05‰BAC.

## ■ Low battery detection

This device using 2 \* 1.5V AAA alkaline batteries, if the batteries are low, “Lob” and the battery symbol in the upper right corner flashes 10 times together with “di-di” sound prompt to remind that need to replace new batteries, then device power off automatically.

Please replace new batteries for normal use.

Low battery prompt

### Warning

**While low battery, device cannot work normally, need to replace new batteries, or it won't power on.**

## ■ Specifications

Vlectro VL-AT

Item	Specifications
Sensor	Semiconductor sensor
Reference Standard	EN16280
Response Time	3-5s
Blowing Time	About 5s
Detection Range	0.00~2.00‰BAC Or 0.000~0.200%BAC Or 0.000 ~1.000mg/L (if ‰BAC and mg/L convection factor is 2.0) i.e. 0.50‰BAC = 0.050%BAC = 0.250mg/L
Sensitivity	0.001% BAC
Sensor Accuracy	+/-0.010 at 0.050% BAC;>0.050% BAC: +/-20%
Optimal Operating Temperature	0°C~50°C
Storage Temperature	-20°C~70°C
Display	About 1.5 inch LCD(positive)
Power Supply	2 * AAA alkaline batteries
Dimensions	110*61*24mm
Net Weight	54g

## ■ Packing accessories

Packing accessories may vary by different countries or operators, please check the packing list according to the order contract.

Standard packing:

1\* VL-AT device

**Vlectro VL-AT**

1\* user manual

1\* color box

Mouthpieces (Quantity optional)

Optional accessories:

2 \* AAA alkaline batteries

**Tip:** If find that the instrument is damaged or missing accessories, please contact your supplier.

## ■ Common faults and solutions

Faults	Possible causes	Solutions
Cannot power on	Low battery.	Replace new batteries.
	Batteries are upside down.	Re-install the batteries in right way.
The two test results differ greatly in a short time	The vent outlet maybe blocked during one test.	Don't block the vent outlet while doing breath test.
	After drinking alcohol, the waiting time is not enough before the alcohol test.	To get a more accurate test result, please abstain from smoking and drinking for at least 20 minutes before testing.
	Before the test, no water or non-alcoholic beverages were used to rinse the mouth and a large amount of alcohol remained in the mouth.	Rinse out your mouth with water or non-alcoholic drinks before using the unit. To ensure the test result more closer to human Blood Alcohol Concentration. Since residues of alcohol in mouth will affect the measured values.

## ■ Maintenance

- ◆ Keep the device (especially the blow hole of the device) away from organic solvents (such as alcohol, etc.) while cleaning the surface of instrument.
- ◆ Don't let alcohol or any other liquid enter into the device, or it will cause permanent damage to the sensor.

Vlectro VL-AT

- ◆ Store the device in right packing box, avoid exposure to sunlight and humidity.

## ■ Precautions

- ◆ During the breath sampling, do NOT block the outlet to avoid a inaccurate test result.
- ◆ Many factors (e.g. environment gas residue, smoke & various smell in human mouth, etc···) may affect the Blood Alcohol Concentration (BAC) level, so the Readings are For Reference Only.
- ◆ Tolerance to alcohol varies from one person to another, and can be affected by such factors as genetics, constitution, gender, percentage of body fat, adaptation to chronic alcohol use, and synergistic effects of drugs. same amount of alcohol, different person the Blood Alcohol Content maybe different, so the test result may be different. As there's a absorption and digestion progress of human's blood alcohol content after drinking, your blood alcohol content (BAC) may change within minutes.

## ■ Important notes

- ◆ The frequency of use and test concentration of device, will affect the sensitivity and life of the sensor. The more frequently used and higher alcohol concentration tested, sensor sensitivity will be reduced and service life will be shortened.
- ◆ Test results for reference only, do not use this device to determine if you can drive or not.
- ◆ The test results may be inadmissible in court.
- ◆ Please bear in mind that this device is designed to measure your blood alcohol concentration/content(BAC),not to estimate how much alcohol you have drunk.
- ◆ For the health and safety of you and others, please don't drive after drink.

## Warranty

- ◆ Vlectro is a registered trademark of Confinity NV. We offer a 2 year warranty on this product.

### Warning

**Neither the re-seller nor manufacturer accepts any liability arising from the use of the unit. Don't drink before you drive!**

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